



## Dispensador para tapa de botella

**Código:** N/D

**Marca:** Belltronic

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

## Detalle

Bottle-Top Dispenser is a device used for liquid separation. This equipment is suitable for general acids and bases, as well as low-concentration strong acids, strong bases and salts.

### Características

1. Excellent chemical resistance.
2. Fully autoclavable at 121°C.
3. Adjustable control knob makes setting the volume easy.
4. Made of PTFE, FEP, BSG, PP.
5. Vapor pressure Max.500mbar, viscositymax. 500mm<sup>2</sup>/s, temperature max.40°C, density max.2.2g/cm<sup>2</sup>.

## Especificaciones técnicas

Compara variantes usando especificaciones técnicas validadas.

| Especificación                | 1 ml - No aplica - $\pm 2,00\%$ , 20 $\mu$ l - modelo 1 | 0,5~5 ml - 0,1 ml - $\pm 0,50\%$ , 25 $\mu$ l - modelo 2 | 1~10 ml - 0,2 ml - $\pm 0,50\%$ , 50 $\mu$ l - modelo 3 | 2,5~25 ml - 0,5 ml - $\pm 0,50\%$ , 125 $\mu$ l - modelo 4 | 5~50 ml - 1,0 ml - $\pm 0,50\%$ , 250 $\mu$ l - modelo 5 |
|-------------------------------|---------------------------------------------------------|----------------------------------------------------------|---------------------------------------------------------|------------------------------------------------------------|----------------------------------------------------------|
| Rango de volumen              | 1 ml                                                    | 0,5~5 ml                                                 | 1~10 ml                                                 | 2,5~25 ml                                                  | 5~50 ml                                                  |
| Subdivisión                   | No aplica                                               | 0,1 ml                                                   | 0,2 ml                                                  | 0,5 ml                                                     | 1,0 ml                                                   |
| Exactitud                     | $\pm 2,00\%$ , 20 $\mu$ l                               | $\pm 0,50\%$ , 25 $\mu$ l                                | $\pm 0,50\%$ , 50 $\mu$ l                               | $\pm 0,50\%$ , 125 $\mu$ l                                 | $\pm 0,50\%$ , 250 $\mu$ l                               |
| Precisión                     | 1,00%, 10 $\mu$ l                                       | 0,10%, 5 $\mu$ l                                         | 0,10%, 10 $\mu$ l                                       | 0,10%, 25 $\mu$ l                                          | 0,10%, 50 $\mu$ l                                        |
| Embalaje                      | 20 uds/cartón                                           | 20 pzs/cartón                                            | 20 unidades/cartón                                      | 20 piezas/cartón                                           | 20 pzas/cartón                                           |
| Peso bruto (kg)               | 14                                                      | 14                                                       | 14                                                      | 14                                                         | 14                                                       |
| Dimensiones del embalaje (mm) | 595*430*275 mm                                          | 595*430*275 mm                                           | 595*430*275 mm                                          | 595*430*275 mm                                             | 595*430*275 mm                                           |

## Aplicaciones

Use la guía para elegir por rango de dosificación y control de precisión. Seleccione cruzando Rango de volumen, Precisión; esos son los ejes que diferencian la cotización técnica.

01

## Rango de volumen

Revise este eje con valores como 1 ml, 0,5~5 ml, 1~10 ml, 2,5~25 ml y 5~50 ml.

02

## Precisión

Revise este eje con valores como 1,00%, 10 µl, 0,10%, 5 µl, 0,10%, 10 µl, 0,10%, 25 µl y 0,10%, 50 µl.

## Cómo elegir rápido

Priorice el volumen que realmente dosifica y después valide graduación, precisión y repetibilidad. Confirme la variante por las especificaciones visibles antes de cotizar.

**Siguiente paso:** [ver producto en Cotecno](#) | [solicitar cotización](#)