


Supplementary Information of Rocker's Accessories

Silicone Stopper Size



Part No.	Item	A Top	B Bottom	C Height	Hole
167110-18	#18 Silicone stopper	Φ 50.5 mm	Φ 44 mm	37 mm	Φ 9.2 mm
167110-17	#16 Silicone stopper	Φ 44 mm	Φ 37.5 mm	35.5 mm	Φ 9.2 mm
167110-16	# 8 Silicone Stopper	Φ 41 mm	Φ 33.4 mm	26 mm	Φ 9.2 mm
167110-16-1	# 8 Silicone Stopper	Φ 41 mm	Φ 33.4 mm	26 mm	Φ 15.6 mm
167110-05	# 5 Silicone Stopper	Φ 29 mm	Φ 24.5 mm	25.5 mm	Φ 9.6 mm

