

VACUUM FILTRATION SYSTEM

ROCKER 300 – MF 31



Features

OIL-FREE & CARE-FREE

Rocker pumps are driven by piston that needs no lubricant, making them maintenance-free and oil-mist-free.

CLAMP-FREE WITH MAGNET

MF 31 bottles are equipped with patented magnetic MF filter holders. Its strong magnetic seal enables one-hand and clamp-free operation.

DRAIN OUT DIRECTLY

MF 31 can drain the filtrate away directly without disassembly.

INTERNATIONAL CERTIFICATION

- CE certification
- UKCA certification
- CSA certification
- IP30 protection
- RoHS certification

Applications

- Suspended solid test
- Vacuum filtration

Model	Rocker 300 – SF 10	
VOLTAGE	110V/60Hz	220V/50Hz
Max. Power	65W	60W
Max. Current	0.7A	0.3A
Max. Vacuum	99 mbar	
Max. Flow Rate	23 L/min	20 L/min
Motor Rotation Speed	1750 RPM	1450 RPM
Horse Power	1/8 HP	
Noise Level	50 dB	
Hose Barb	ID8 (5/16 inch)	
Vacuum Regulator	Yes	
Overflow Protection	Yes	
Thermal Protection	Yes	