

ROCKER SOLVENT PURIFICATION SYSTEM



**ROCKER 300C
VF 12**

Features

OIL-FREE & CARE-FREE

Rocker 300C / Chemker 300 / Lafil 300C pumps are driven by diaphragm that needs no lubricant, making them maintenance-free and oil-mist-free.

CHEMICAL RESISTANT

Rocker 300C / Chemker 300 / Lafil 300C - VF 12 are resistant in different degrees to wide range of chemicals and solvents.

ONE-STEP FILTRATION

Filter and store directly in all GL45 bottles to prevent repetitive transfers and cross-contamination

INTERNATIONAL CERTIFICATION

- CE certification
- IP 30 protection (Rocker 300C)
- RoHS certification (Rocker 300C and Chemker 300)

Applications

- Solvent purification
- HPLC mobile phase purification
- Vacuum filtration



**CHEMKER 300
VF 12**



**LAFIL 300C
VF 12**

Model	Rocker 300C - VF 12		Chemker 300 - VF 12		Lal 300C - VF 12	
VOLTAGE	110V	220V	110V	220V	110V	220V
Max. Power	65W	60W	70W	60W	65W	60W
Max. Vacuum	99 mbar abs.		99 mbar abs.		99 mbar abs.	
Max. Flow Rate	23 L/min @60Hz	20 L/min @50Hz	20 L/min @60Hz	18 L/min @50Hz	23 L/min @60Hz	20 L/min @50Hz
Vacuum Regulator	Yes		-		Yes	
Noise Level	50 dB		50 dB		50 dB	
Model of Apparatus	VF 12					
Funnel	300 mL					
Storage Bottle	1000 mL					
Filter Diameter	47 mm					