



Extractor medico de plasma sanguineo para bolsas de 100 a 400 ml

Código: B2-01-010591

Marca: Belltronic

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Plasma Extractor is designed to extract blood component from centrifuged bags. It is widely used in hospital, blood station, etc. After the plasma of blood bags is layered, the pressure of plasma-separating clip enables different components of blood to be extruded from top to bottom.

Features

- 1. Manual system
- 2. Easy to observe the separation of the plasma and red blood cells

Especificaciones técnicas

Especificación	20N - 100ml, 200ml, 300ml, 400ml - Acero inoxidable
Fuerza de sujeción	20N
Bolsa de sangre aplicable	100ml, 200ml, 300ml, 400ml
Material	Acero inoxidable
Dimensiones externas (An*Pr*Al)	165*250*250 mm
Dimensiones del embalaje (An*Pr*Al)	310*310*240mm
Peso neto	2.5kg
Peso bruto	2.8kg

Aplicaciones

1. Confirm your bags

Confirm that the volume of your blood bags is among the applicable ones.

- **Applicable Blood Bag:** 100ml, 200ml, 300ml, 400ml

2. Plan the space

Plan the bench space for the extractor.

- **External Size (W*D*H):** 165*250*250 mm



Importadora Gyb Limitada
Padre orellana 1476
Santiago CL13101 Chile

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)