



Microscopio biológico ocular P10×/18 mm objetivo 4X AN = 0,10 WD = 37,50 mm

Código: B2-01-010649

Marca: Belltronic

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Introducción

El microscopio biológico B2-01-010649 es su puerta de entrada para explorar los detalles más minúsculos de la vida con precisión y claridad inigualables. Este microscopio de vanguardia está diseñado para satisfacer las necesidades tanto de principiantes como de investigadores experimentados, lo que lo convierte en una herramienta esencial para diversas disciplinas científicas y entornos educativos.

Especificaciones técnicas

Compara variantes usando especificaciones técnicas validadas.

Especificación	P10X/18mm, P16X/11mm - 4X N.A=0.10 WD=37.50mm - Yellow, green, blue
View Head	Compensation free binocular head, Inclined at 30°, 360° rotatable, Interpupillary 55~75mm
Eyepiece	P10X/18mm, P16X/11mm
Achromatic Objective	4X N.A=0.10 WD=37.50mm
Focusing	Coaxial Coarse/Fine Adjustable System, Range 15mm, Fine Division 0.002mm
Stage	Mechanical Stage 130*130mm, Moving Range 70*30mm, Vernier 0.1mm
Condenser	Abbe NA1.25 condenser with Iris Diaphragm/Filter
Illumination	LED dual-color light source, adjustable brightness
Filter	Yellow, green, blue
Power Supply	AC85~240V, 50/60Hz
Package Size (W*D*H)	435*330*235mm
Gross Weight	5.5kg

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)