



Maquina de ensayo de compresion 2000 kN, espacio 320 mm

Código: B2-01-011018

Marca: Belltronic

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Digital Display Compression Testing Machine mainly used for concrete and other material compression resist testing Widely used in metallurgy, building materials, space flight and aviation, colleges and universities, R&D institution lines. The test operation and data processing meets standard requirements

Especificaciones técnicas

Especificación	4% 100%FS(80KN 2000KN) - 320mm - 250*250mm(La*An)
Carga máxima	2000KN
Clase	Clase 1
Ajuste del espacio de ensayo	Ajuste manual por husillo de bolas
Modo de control	Control manual de carga
Estructura del gabinete de control de la fuente de aceite	Gabinete de control de fuente de aceite de mesa en una pieza
Rango de fuerza de ensayo	4% 100%FS(80KN 2000KN)
Exactitud del valor indicado	±1%
Carrera del pistón	50mm
Separación máxima entre placas de compresión superior e inferior	320mm
Separación máxima entre columnas	270mm
Tamaño de las placas de compresión	250*250mm(La*An)
Pantalla	Pantalla digital LCD
Dimensiones del equipo principal	900 400*1150mm
Alimentación eléctrica	1KW I AC 380V, 50Hz(se puede personalizar)
Peso	750kg

Aplicaciones

1. Confirm your specimens

Confirm that your specimens fit between the compression plates and between the columns.

- **The Largest Spacing Of Upper And Lower Pressure Plate:** 320mm
- **The Largest Spacing Of Two Columns:** 270mm

2. Review the accuracy

Review the indicated value accuracy against your test standard.

- **Indicating Value Accuracy:** $\pm 1\%$

3. Plan the installation

Plan the space and floor strength for the weight of the machine.

- **The Host Dimension:** 900 400*1150mm
- **Weight:** 750kg

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)