



Armario de secado médico

Código: B2-01-011010

Marca: Belltronic

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Medical drying cabinet is a large-capacity medical drying equipment developed by our company according to the needs of medical and scientific research institutes. It is used to dry reusable medical products such as surgical instruments, anesthesia instruments, glassware, and catheter endoscopes. The operation is simple and flexible, the drying temperature setting range is 40°C~80°C, and the drying time setting range is 1min~999min. Users can freely set drying procedures suitable for medical institutions, scientific research institutes, laboratories and other medical, scientific research and production units according to their needs.

Características

1. Made of high-quality 304 stainless steel, with low bacterial adhesion rate, safe and convenient to use.
2. Insulation interlayer, effectively prevent heat loss.
3. With a glass window, you can observe the internal working conditions at any time.
4. Reasonable air duct and circulation system to ensure good temperature uniformity.
5. Shelf height is adjustable.
6. Using brand motor and fan blades, the wind speed is strong, which is conducive to the renewal and circulation of air in the inner cavity.
7. It can centrally handle a large number of required drying items, and can simultaneously dry surgical instruments, anesthesia pipes, and respiratory treatment items.
8. inch PLC human-computer interaction control system, flexible and convenient operation, safe and stable work.

Especificaciones técnicas

Compara variantes usando especificaciones técnicas validadas.

Especificación	400 - 1~999min - 2*soportes para frascos de humidificación - modelo 1	600 - 1~999 min - 2*soportes para botellas de humidificación - modelo 2
Capacidad (L)	400	600
Material	S30408	S30408
Puerta	Puerta visual hermética de vidrio templado de doble apertura	Puerta visual hermética de vidrio templado de doble apertura
Rango de temperatura	40~80°C	40~80°C
Ajuste de tiempo	1~999min	1~999 min
Alimentación eléctrica	CA 220 V, 50/60 Hz (estándar); CA 380 V, 50/60 Hz (trifásico, opcional)	AC380V, 50/60Hz (estándar); AC220V, 50/60Hz (monofásico, opcional)
Consumo	5,8 kW	6,6 kW
Accesorios estándar	9*estantes; 9*bandejas; 2*soportes para catéter	9*estantes; 9*bandejas; 2*soportes para catéteres

Especificación	400 - 1~999min - 2*soportes para frascos de humidificación - modelo 1	600 - 1~999 min - 2*soportes para botellas de humidificación - modelo 2
Accesorios opcionales	2*soportes para frascos de humidificación	2*soportes para botellas de humidificación
Dimensiones externas (A×P×H)	950*550*2170mm	950*750*2170mm
Dimensiones del embalaje (A×P×H)	1060*870*2350mm	1060*1050*2350mm
Peso neto	320kg	350 kg
Peso bruto	380kg	420 kg

Aplicaciones

Armario de secado médico: elija la variante comparando Capacidad (L).

Que cambia entre variantes

B2-01-011010.1: Capacidad (L): 400.

B2-01-011010.2: Capacidad (L): 600.

Escenarios de uso recomendados

Use esta alternativa si necesita capacidad (l): 400.

Use esta alternativa si necesita capacidad (l): 600.

Cómo elegir rápido

Revise primero Capacidad (L); si esos puntos calzan con el proceso, confirme los accesorios y condiciones de instalacion.

Recomendación de compra

Seleccione la alternativa que cubra Capacidad (L).

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)