

Orbital Shaker

Features

Orbital Shaker one powerful variable speed shaker which provides efficient orbital motion. Speed and time are under microprocessor control. 250-12000 RPM shaker can be applicable in different lab: microbiology, chemistry, immunology, biochemistry and molecular biology







Specification

Model	BellPlaUP-100	BellPlaUP-1000
Shaking Speed	50-1200rpm	50-1200rpm
Orbit	10mm	20mm
Timer range	1min-9999min	1min-9999min
Load capacity	See below table	See below table
Max. capacity	2.0kg	2.0kg

Plataforma Universal, BellPlaUP-1

Código: B-01-29-000402-0105

Marca: Belltronic

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Plataforma Universal. Ficha generada desde el nombre comercial.

1. Confirmá el nombre: Plataforma Universal, BellPlaUP-1.
2. Validá con ficha de proveedor si necesitás norma.

Tipo: Plataforma Universal

Referencia de nombre: Plataforma Universal, BellPlaUP-1

Dato del nombre: 1

Origen de ficha: NAME_ONLY (datos del nombre comercial)

Uso previsto: Laboratorio / ensayo / control de calidad

Especificaciones técnicas

Aplicaciones

Uso de laboratorio / ensayo

Verificar compatibilidad con protocolo

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)