



Sistema Minneapolis Blower Door™ con Manómetro DG-1000

Código: BD3-KIT-001

Marca: The Energy Conservatory

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Detalle técnico

Description Minneapolis Blower Door TM System with DG-1000 The Minneapolis Blower Door TM has long been recognized as the best designed and supported building airtightness testing system in the world. Combined with specialized accessories and testing procedures developed by The Energy Conservatory, the Minneapolis Blower Door TM is the system of choice for utility DSM programs, energy raters, HVAC contractors and weatherization professionals. Please contact TEC to place your order if you need custom fan calibration service for specific testing standards. Minneapolis Blower Door TM System Features Includes everything required for fully automated testing Durable injection molded fan housing with a 5 year TEC repair warranty Adjustable aluminum frame allows for quick and easy set-up Versatile accessory case features padded pockets for laptops and files Provides accurate flow measurements...

Marca: The Energy Conservatory. **Proveedor:** THE ENERGY CONSERVATORY INC..

Especificaciones técnicas

Compara modelos principales con especificaciones tecnicas validadas.

Especificación	Sistema o kit de ensayo - Manómetros Digitales
Tipo de producto	Sistema o kit de ensayo
Subfamilia Cotecno	Manómetros Digitales
Línea TEC	Minneapolis Blower Door™ / Kits Minneapolis Blower Door™
Peso declarado por proveedor	1 lbs
Documentación disponible	Formulario de servicio de equipo; Folleto general; Documento de referencia; Folleto técnico;

Aplicaciones

- Medición de presión
- Lectura de caudal
- Registro de datos de ensayo

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)