

**ACCESSORIES**

For Molds, please see page 6.7.

Length Comparator, Dial Indicator

ASTM C151, C157, C227, C311, C341, C342, C452, C490, C596; AASHTO M210, T107, T160.

Length comparators measure length changes of hardened cement paste, mortar and concrete prismatic specimens. Indicators are mounted on a sturdy upright support attached to a solid triangular base. All units include stationary and movable anvils designed to fit H-3260 gauge studs, which are cast into test specimens and an invar reference bar. Comparators for 10" (254mm) or 16" (406mm) effective length samples with 3.625" (92mm) dia. dial gauge with a range of .400" and gradations of .0001". The Dial is marked 0-10. Includes stationary and movable anvils designed to fit the H-3260 gauge studs that are cast on test specimens and an invar reference bar. Unit accommodates test specimens up to 4" x 4" (102 x 102mm) cross section.

Length Comparator, 10" Effective Length H-3250**Length Comparator, 16" Effective Length** H-3248

Ship wt. 32.7lbs. (14.8kg)

Length Comparator, Digital Indicator

ASTM C151, C157, C227, C311, C341, C342, C452, C490, C596; AASHTO M210, T107, T160.

Length comparators measure length changes of hardened cement paste, mortar and concrete prismatic specimens. Indicators are mounted on a sturdy upright support attached to a solid triangular base. All units include stationary and movable anvils designed to fit H-3260 gauge studs, which are cast into test specimens and an invar reference bar. Comparators for 10" (254mm) or 16" (406mm) effective length samples with digital indicator with a range of .600" and resolution of .0001". Measures inches and millimeters includes batteries. Unit can be zeroed at any point on the range and can be switched from inches to mm by pressing a button. Includes stationary and movable anvils designed to fit the H-3260 gauge studs that are cast on test specimens and an invar reference bar. Unit accommodates test specimens up to 4" x 4" (102 x 102mm) cross section.

Length Comparator, 10" Length H-3250D**Length Comparator, 16" Length** H-3248D

Ship wt. 36.5lbs. (16.8kg)

Comparator Replacement Components

Description	Model
Invar reference bar for use with H-3250 and 10" (254mm) specimens w/ sleeve	H-3249A
Invar reference bar for use with H-3248 and 16" (406mm) specimens w/ sleeve	H-3249A.16
Invar Sleeve	H-3249.2
Dial Indicator w/Collar & Anvil	H-3250.3
Digital Indicator, Battery Powered	H-3250.3DA
Anvil, Top with collar (one each)	H-3250.4
Anvil, Top only	H-3250.4.4
Collar w/ set screw	H-3250.4.3
Modified Collar with Set Screw	H-3250.4.6
Elevating screw, nut collar & anvil assembly, Bottom	H-3250.7
Adapter for 5" specimens	H-3250.8

Custom size Invar Reference Bars are available, please inquire: 1.800.544.7220.

NOTES

Length Comparators can be custom made for any length prism molds.