




Sterlitech Sterile Cellulose Nitrate, MCE Membranas

Código: N/D

Marca: STERLITECH

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Sterlitech Sterile Cellulose Nitrate, MCE Membranas

Its features include high porosity, high protein binding can be blocked by pre-treatment or utilized in application, Triton-free, biologically inert with good thermal stability, and high degree of internal surface area for greater adsorption of product.

Especificaciones técnicas

Todos bajo pedido por ahora. Filtra por poro, tamaño o pack para elegir el SKU.

Especificación	0.22 µm · 47 mm · pack 100 · MFCN047022GS	0.22 µm · 47 mm · pack 100 · MFCN047022S	0.45 µm · 47 mm · pack 100 · MFCN047045GWS	0.45 µm · 47 mm · pack 100 · MFCN047045S	0.8 µm · 47 mm · pack 100 · MFCN047080GWS	0.8 µm · 47 mm · pack 100 · MFCN047080WS	1 µm · 47 mm · pack 100 · MFCN047100GWS	1 µm · 47 mm · pack 100 · MFCN047100S
Tamaño de poro	0.22	0.22	0.45	0.45	0.8	0.8	1	1
Cuadrícula	Yes	No	Yes	No	Yes	No	Yes	No
Diámetro (mm)	47	47	47	47	47	47	47	47
Rejilla	Yes	No	Yes	No	Yes	No	Yes	No
Tamaño de paquete	100	100	100	100	100	100	100	100

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)