



Kit de O-ring Micro Kleansep EPDM

Código: 1201672

Marca: STERLITECH

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle Ceramic Housings

Producto publicado con base en la ficha vigente del proveedor Sterlitech.

SKU	Nombre proveedor	Fabricante	Categoría proveedor	Fuente
1201672	O-ring/Gasket Kit for Micro Kleansep, EPDM	Micro-Kleansep™	Filters /// Membrane Disc Filters Membrane/Process Development /// Membrane Housings and Elements	Ficha proveedor

Especificaciones técnicas

Familia Odoo	Kit de O-ring Micro Kleansep EPDM
Fichas fuente	1

1201672

SKU Sterlitech	1201672
Ficha del fabricante	O-ring/Gasket Kit for Micro Kleansep, EPDM
URL fuente	

Datos visibles del fabricante

O-ring/Gasket Kit for Micro Kleansep, EPDM SKU: 1201672 Brand: Micro-Kleansep™

Contenido visible del fabricante

Applications Applications Concentration of biological media for the recovery of antibodies, enzymes, proteins, etc. Water rendered apyrogenic Diafiltration Separation or clarification of molecules of very distinct sizes Cell recovery and washing for the concentration of viruses, bacteria, yeasts and animal cells Purification of protein solutions Elimination of cell fermentation and culture metabolites by continuous fermentation to improve yields and optimise culture cycles Dialysis/desalination of fragile proteins by low-pressure operations (< 1 bar) Organic Solvent Nanofiltration (OSN) Wide range of feasibility studies in chemistry, waste treatment, etc.

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)