

	CF016	CF042	Seipa
Active Membrane Area	20.8 cm ² (3.2 in ²)	42 cm ² (6.5 in ²)	140 cm ² (22 in ²)
Hold-Up Vol.	13 mL (0.44 oz)	17 mL (0.57 oz)	70 mL (2.4 oz)
Connections	Feed: $\frac{1}{8}$ in FNPT (base of cell) Concentrate: $\frac{1}{8}$ in FNPT (base of cell) Filtrate: $\frac{1}{8}$ in FNPT (top of cell)	Feed: $\frac{1}{8}$ in FNPT (base of cell) Concentrate: $\frac{1}{8}$ in FNPT (base of cell) Filtrate: $\frac{1}{8}$ in FNPT (top of cell)	Feed: $\frac{1}{8}$ in FNPT (base of cell) Concentrate: $\frac{1}{8}$ in FNPT (base of cell) Filtrate: $\frac{1}{8}$ in FNPT (top of cell)
Outer Dimensions	12.7 x 10 x 8.3 cm (5 x 4 x 3.25 in)	12.7 x 10 x 8.3 cm (5 x 4 x 3.25 in)	16.51 x 21.3 x 5 cm (6.5 x 8.38 x 2.01 in)
Active Area Dimensions	4.52 x 4.52 cm (1.78 x 1.78 in)	9.207 x 4.572 cm (3.625 x 1.8 in)	9.7 x 14.7 cm (3.81 x 5.78 in)
Slot Depth	0.23 cm (0.09 in)	0.23 cm (0.09 in)	0.19 cm (0.075 in)
Slot Width	39 mm (1.54 in)	39 mm (1.54 in)	95.3 mm (3.75 in)

Kits de celdas de flujo cruzado y tangencial

Código: N/D

Marca: STERLITECH

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Kits de celdas de flujo cruzado y tangencial

Kits y accesorios para ensayos de flujo cruzado/tangencial con membranas planas. Se utilizan para evaluar condiciones de proceso y comportamiento de permeado y retentado.

Especificaciones técnicas

Todos bajo pedido por ahora. Filtra por poro, tamaño o pack para elegir el SKU.

Especificación	1201805	1202012	1202013	1240312
Tipo	[1201805] Kits de celdas de flujo cruzado y tangencial (1201805)	[1202012] Kits de celdas de flujo cruzado y tangencial (1202012)	[1202013] Kits de celdas de flujo cruzado y tangencial (1202013)	[1240312] Kits de celdas de flujo cruzado y tangencial (1240312)

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)