



## Laboratorio Pumps

**Código:** N/D

**Marca:** STERLITECH

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

## Detalle Laboratory Pumps

Producto publicado con base en la ficha vigente del proveedor Sterlitech.

| SKU       | Nombre proveedor                                                                                | Fabricante | Categoría proveedor                       | Fuente          |
|-----------|-------------------------------------------------------------------------------------------------|------------|-------------------------------------------|-----------------|
| 1220173   | PILOT Service Kit - Diaphragm Replacement Parts (TLD3000)                                       | Sterlitech | Laboratory Equipment /// Laboratory Pumps | Ficha proveedor |
| 1220174   | PILOT Service Kit - Major Rebuild Parts (TLD3000)                                               | Sterlitech | Laboratory Equipment /// Laboratory Pumps | Ficha proveedor |
| 1220175   | PILOT Service Kit - Diaphragm Replacement Parts (TLD5000)                                       | Sterlitech | Laboratory Equipment /// Laboratory Pumps | Ficha proveedor |
| 1220176   | PILOT Service Kit - Major Rebuild Parts (TLD5000)                                               | Sterlitech | Laboratory Equipment /// Laboratory Pumps | Ficha proveedor |
| 1220177   | PILOT Service Kit - Rubber Suction Cup Feet (4/pk) (TLD3000/5000)                               | Sterlitech | Laboratory Equipment /// Laboratory Pumps | Ficha proveedor |
| 1240079   | PILOT Pressure Regulator Kit Assembly                                                           | Sterlitech | Laboratory Equipment /// Laboratory Pumps | Ficha proveedor |
| 167801-07 | Filter Cartridge for Rocker 600/610/801/811 - without Prote and overflow protection (167801-07) | Rocker     | Laboratory Equipment /// Laboratory Pumps | Ficha proveedor |

## Especificaciones técnicas

Todos bajo pedido por ahora. Filtra por poro, tamaño o pack para elegir el SKU.

| Especificación | 1220173                               | 1220174                               | 1220175                               | 1220176                               | 1220177                               | 1240079                               | 167801-07                                 |
|----------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|-------------------------------------------|
| Tipo           | [1220173] Laboratorio Pumps (1220173) | [1220174] Laboratorio Pumps (1220174) | [1220175] Laboratorio Pumps (1220175) | [1220176] Laboratorio Pumps (1220176) | [1220177] Laboratorio Pumps (1220177) | [1240079] Laboratorio Pumps (1240079) | [167801-07] Laboratorio Pumps (167801-07) |

**Siguiente paso:** [ver producto en Cotecno](#) | [solicitar cotización](#)