

Rotor angular 48 x 1,5/2 ml para centrifuga de mesa refrigerada de alta velocidad

Código: B2-01-010384-R05

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Angle Rotor

Especificaciones técnicas

Especificación	-R05
Tipo de rotor	Rotor angular
Volumen	48*1.5/2ml
Velocidad máxima (rpm)	13000
Fuerza centrífuga relativa máxima (xg)	20122
Adaptadores	0.2ml, 0.5ml

Aplicaciones

1. Review your tubes

Compare the number and volume of the tubes you process per run with the rotor capacity.

- **Volume:** 48*1.5/2ml

If you work with smaller tubes, confirm that the listed adapters match them.

- **Adapters:** 0.2ml, 0.5ml

2. Review speed and centrifugal force

Compare the speed and relative centrifugal force your protocol requires with the rotor maximums.

- **Max. Speed(rpm):** 13000
- **Max. RCF(xg):** 20122

3. Confirm your centrifuge

Confirm that your centrifuge accepts this rotor type before buying it.

- **Rotor Type:** Angle Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)