

Rotor angular 12 x 10/15 ml de baja velocidad para centrifuga de mesa refrigerada de alta velocidad

Código: B2-01-010384-R12

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Angle Rotor

Especificaciones técnicas

Especificación	-R12
Tipo de rotor	Rotor angular
Volumen	12*10/15ml(baja velocidad)
Velocidad máxima (rpm)	6000
Fuerza centrífuga relativa máxima (xg)	4387
Adaptadores	Tubo de recolección de sangre de 2/3/4/5/7ml, tubo de recolección de sangre de 5/7ml, 15ml(▲)

Aplicaciones

1. Review your tubes

Compare the number and volume of the tubes you process per run with the rotor capacity.

- **Volume:** 12*10/15ml(low speed)

If you work with blood collection tubes, confirm that the listed adapters match them.

- **Adapters:** 2/3/4/5/7ml blood collection tube, 5/7ml blood collection tube, 15ml(▲)

2. Review speed and centrifugal force

This is a low-speed rotor. Compare the speed and relative centrifugal force your protocol requires with the rotor maximums.

- **Max. Speed(rpm):** 6000
- **Max. RCF(xg):** 4387

3. Confirm your centrifuge

Confirm that your centrifuge accepts this rotor type before buying it.

- **Rotor Type:** Angle Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)