

Rotor angular 4 x 100 ml para centrifuga de mesa refrigerada de alta velocidad

Código: B2-01-010384-R18

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Angle Rotor

Especificaciones técnicas

Especificación	-R18
Tipo de rotor	Rotor angular
Volumen	4*100ml(▲/●)
Velocidad máxima (rpm)	12000
Fuerza centrífuga relativa máxima (xg)	17700
Adaptadores	—

Aplicaciones

1. Confirm the tubes you will spin

Compare the rotor's number of positions and tube volume with the samples you process in each run.

- **Volume:** 4*100ml(▲/●)

2. Check speed and centrifugal force

Check that the rotor's maximum speed and maximum relative centrifugal force meet what your separation protocol requires.

- **Max. Speed(rpm):** 12000
- **Max. RCF(xg):** 17700

3. Choose the rotor type

Choose between an angle rotor and a swing rotor according to the separation you need, and confirm that the rotor matches your centrifuge.

- **Rotor Type:** Angle Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)