

Rotor oscilante 8 x 100 ml para centrifuga de mesa refrigerada de alta velocidad

Código: B2-01-010384-R32

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Swing Rotor

Especificaciones técnicas

Especificación	-R32
Tipo de rotor	Rotor oscilante
Volumen	8*100ml
Velocidad máxima (rpm)	4000
Fuerza centrífuga relativa máxima (xg)	3077
Adaptadores	32*2/3/4ml (tubo de recolección de sangre), 32*5/7ml (tubo de recolección de sangre)

Aplicaciones

1. Confirm the capacity and the tubes

Compare the number and volume of tubes you process per run with the volume of the rotor.

- **Volume:** 8*100ml

Check whether the listed adapters accept the tubes used in your laboratory.

- **Adapters:** 32*2/3/4ml blood collection tube, 32*5/7ml blood collection tube

2. Check the speed and compatibility

Check that the maximum speed and relative centrifugal force of the rotor cover what your separation protocol requires.

- **Max. Speed(rpm):** 4000
- **Max. RCF(xg):** 3077

Confirm that the rotor type matches the centrifuge you have and your separation method.

- **Rotor Type:** Swing Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)