

Rotor oscilante 4 x 250 ml de bioseguridad para centrifuga de mesa refrigerada de alta velocidad

Código: B2-01-010384-R33

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Swing Rotor

Especificaciones técnicas

Especificación	-R33
Tipo de rotor	Rotor oscilante
Volumen	4*250ml(bioseguridad)
Velocidad máxima (rpm)	4000
Fuerza centrífuga relativa máxima (xg)	3578
Adaptadores	48*1.5/2ml, 40*2/3/4/5/7ml (tubo de recolección de sangre), 48*2/3/4/5/7ml (tubo de recolección de sangre), 24*10ml(▲), 28*10ml, 20*15ml, 4*50ml(▲), 8*50ml(●), 4*100ml

Aplicaciones

1. Confirm the capacity and the tubes

Compare the number and volume of tubes you process per run with the volume of the rotor.

- **Volume:** 4*250ml(biosafety)

Check whether the listed adapters accept the tubes used in your laboratory.

- **Adapters:** 48*1.5/2ml, 40*2/3/4/5/7ml blood collection tube, 48*2/3/4/5/7ml blood collection tube, 24*10ml(▲), 28*10ml, 20*15ml, 4*50ml(▲), 8*50ml(●), 4*100ml

2. Check the speed and compatibility

Check that the maximum speed and relative centrifugal force of the rotor cover what your separation protocol requires.

- **Max. Speed(rpm):** 4000
- **Max. RCF(xg):** 3578

Confirm that the rotor type matches the centrifuge you have and your separation method.

- **Rotor Type:** Swing Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)