

Rotor angular 24 x 10 ml para centrifuga de mesa refrigerada de alta velocidad

Código: B2-01-010384-R46

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Angle Rotor

Especificaciones técnicas

Especificación	-R46
Tipo de rotor	Rotor angular
Volumen	24*10ml
Velocidad máxima (rpm)	5000
Fuerza centrífuga relativa máxima (xg)	3354
Adaptadores	1.5/2ml, tubo de recolección de sangre de 2/3/4ml

Aplicaciones

1. Capacity and tube type

Confirm the rotor type and that the number and volume of the tubes you process per run match those of the rotor.

- **Rotor Type:** Angle Rotor
- **Volume:** 24*10ml

2. Speed and centrifugal force

Compare the speed and relative centrifugal force your protocol requires with the rotor maximums.

- **Max. Speed(rpm):** 5000
- **Max. RCF(xg):** 3354

3. Adapters for smaller tubes

If you use smaller-volume tubes, check that they appear among the adapters listed for the rotor.

- **Adapters:** 1.5/2ml, 2/3/4ml blood collection tube

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)