

Rotor angular 24 x 10/15 ml para centrifuga de mesa refrigerada de alta velocidad

Código: B2-01-010384-R49

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Angle Rotor

Especificaciones técnicas

Especificación	-R49
Tipo de rotor	Rotor angular
Volumen	24*10/15ml
Velocidad máxima (rpm)	6000
Fuerza centrífuga relativa máxima (xg)	5031
Adaptadores	tubo de recolección de sangre de 2/3/4ml, tubo de recolección de sangre de 5/7ml, 15ml(▲)

Aplicaciones

1. Capacity and tube type

Confirm the rotor type and that the number and volume of the tubes you process per run match those of the rotor.

- **Rotor Type:** Angle Rotor
- **Volume:** 24*10/15ml

2. Speed and centrifugal force

Compare the speed and relative centrifugal force your protocol requires with the rotor maximums.

- **Max. Speed(rpm):** 6000
- **Max. RCF(xg):** 5031

3. Adapters for blood collection tubes

If you use blood collection tubes or other smaller tubes, check that they appear among the adapters listed for the rotor.

- **Adapters:** 2/3/4ml blood collection tube, 5/7ml blood collection tube , 15ml(▲)

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)