

Rotor angular 6 x 50 ml para centrifuga de mesa refrigerada de alta velocidad, 15000 rpm, tubo de 106 mm

Código: B2-01-010440-R07

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Angle Rotor

Especificaciones técnicas

Especificación	-R07
Rotor	Rotor angular
Velocidad máxima (rpm)	15000
Fuerza centrífuga relativa (xg)	23640
Capacidad	6*50ml(●), tubo de fondo redondo
Dimensiones del tubo (mm)	Φ29*106

Aplicaciones

1. Confirm the tubes you will centrifuge

Compare the rotor capacity and the tube dimensions with the tubes used in your laboratory.

- **Capacity:** 6*50ml(●)
- **Dimension of Tube (mm):** Φ29*106

Use round-bottom tubes: the (●) mark next to the capacity indicates that bottom type.

- **Capacity:** 6*50ml(●)

2. Check the speed and force you need

Check that the maximum speed and the relative centrifugal force meet what your separation protocol requires.

- **Max. Speed(rpm):** 15000
- **RCF (xg):** 23640

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)