

Rotor oscilante 72 tubos de recoleccion de sangre de 2/3/4/5/7 ml para centrifuga de piso de baja velocidad

Código: B2-01-012213-R03

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Swing Rotor

Especificaciones técnicas

Especificación	-R03
Tipo de rotor	Rotor oscilante
Volumen	72*2/3/4/5/7ml tubo de recolección de sangre
Velocidad máxima (rpm)	4000
Fuerza centrífuga relativa máxima (xg)	2808
Adaptadores	—

Aplicaciones

1. Tube capacity and type

Compare the number and size of tubes you process per run with the rotor capacity.

- **Volume:** 72*2/3/4/5/7ml blood collection tube

2. Speed and centrifugal force

Check that the speed and relative centrifugal force required by your protocol do not exceed the rotor maximums.

- **Max. Speed(rpm):** 4000
- **Max. RCF(xg):** 2808

3. Centrifuge compatibility

Confirm that your floor-type low-speed centrifuge accepts this swing rotor.

- **Rotors Type:** Swing Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)