

Rotor oscilante 16 x 10/15 ml para centrifuga de piso de baja velocidad

Código: B2-01-012213-R06

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Swing Rotor

Especificaciones técnicas

Especificación	-R06
Tipo de rotor	Rotor oscilante
Volumen	16*10/15ml
Velocidad máxima (rpm)	4000
Fuerza centrífuga relativa máxima (xg)	2957
Adaptadores	2/3/4ml tubo de recolección de sangre, 5/7ml tubo de recolección de sangre, 15ml(▲)

Aplicaciones

1. Tube capacity and type

Compare the number and size of tubes you process per run with the rotor capacity.

- **Volume:** 16*10/15ml

If you will use tubes of another size, check the adapters the rotor accepts.

- **Adapters:** 2/3/4ml blood collection tube, 5/7ml blood collection tube, 15ml(▲)

2. Speed and centrifugal force

Check that the speed and relative centrifugal force required by your protocol do not exceed the rotor maximums.

- **Max. Speed(rpm):** 4000
- **Max. RCF(xg):** 2957

3. Centrifuge compatibility

Confirm that your floor-type low-speed centrifuge accepts this swing rotor.

- **Rotors Type:** Swing Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)