

Rotor oscilante 4 x 100 ml para centrifuga de piso de baja velocidad

Código: B2-01-012213-R14

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Swing Rotor

Especificaciones técnicas

Especificación	-R14
Tipo de rotor	Rotor oscilante
Volumen	4*100ml
Velocidad máxima (rpm)	4500
Fuerza centrífuga relativa máxima (xg)	3962
Adaptadores	16*2/3/4ml tubo de recolección de sangre, 16*5/7ml tubo de recolección de sangre

Aplicaciones

1. Confirm the rotor capacity

Compare the number of tubes and the volume you process per run with the rotor capacity.

- **Volume:** 4*100ml

Check that the rotor type matches the centrifuge you own or plan to acquire.

- **Rotors Type:** Swing Rotor

2. Review the speed and centrifugal force

Confirm that the maximum speed and relative centrifugal force of the rotor cover what your separation protocol requires.

- **Max. Speed(rpm):** 4500
- **Max. RCF(xg):** 3962

3. Choose the adapters for your tubes

Review which blood collection tube adapters the rotor accepts and order the ones that match the tubes you use.

- **Adapters:** 16*2/3/4ml blood collection tube, 16*5/7ml blood collection tube

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)