

Rotor oscilante para centrifuga 4 x 250 ml (bioseguridad)

Código: B2-01-012213-R18

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Swing Rotor

Especificaciones técnicas

Especificación	-R18
Tipo de rotor	Rotor oscilante
Capacidad	4*250ml(bioseguridad)
Velocidad máxima (rpm)	4000
FCR máxima (×g)	3041
Adaptadores	48*1.5/2ml, 40*2/3/4/5/7ml tubo de recolección de sangre, 48*2/3/4/5/7ml tubo de recolección de sangre, 24*10ml(▲), 28*10ml, 20*15ml, 4*50ml(▲), 8*50ml(●), 4*100ml

Aplicaciones

1. Check capacity

Check the stated rotor type and capacity to select the component you need.

- **Rotors Type:** Swing Rotor
- **Volume:** 4*250ml(biosafety)

2. Check performance

Compare the stated maximum speed and maximum RCF with your procedure requirements.

- **Max. Speed(rpm):** 4000
- **Max. RCF(×g):** 3041

3. Check adapters

Check the listed adapters against the tubes you use.

- **Adapters:** 48*1.5/2ml, 40*2/3/4/5/7ml blood collection tube, 48*2/3/4/5/7ml blood collection tube, 24*10ml(▲), 28*10ml, 20*15ml, 4*50ml(▲), 8*50ml(●), 4*100ml

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)