

Rotor oscilante 4 x 750 ml de bioseguridad para centrifuga de piso de baja velocidad

Código: B2-01-012213-R24

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Swing Rotor

Especificaciones técnicas

Especificación	-R24
Tipo de rotor	Rotor oscilante
Volumen	4*750ml(bioseguridad)
Velocidad máxima (rpm)	4000
Fuerza centrífuga relativa máxima (xg)	3399
Adaptadores	108*2/3/4/5/7ml tubo de extracción de sangre, 60*15ml(●), 60*15ml(▲), 24*50ml(▲), 24*50ml(●), 28*50ml(●), 28*50ml(▲) (●: tubo de fondo redondo; ▲: tubo de fondo cónico)

Aplicaciones

1. Capacity and tubes

Confirm how many samples you process per run and the volume of each container, and compare it with the rotor capacity.

- **Volume:** 4*750ml(biosafety)

Check whether your blood collection tubes or other tubes need any of the listed adapters.

- **Adapters:** 108*2/3/4/5/7ml blood collection tube, 60*15ml(●), 60*15ml(▲), 24*50ml(▲), 24*50ml(●), 28*50ml(●), 28*50ml(▲)

In the adapters, (●) marks a round-bottom tube and (▲) marks a sharp-bottom (conical) tube; choose the adapter that matches the bottom of your tubes.

- **Adapters:** 108*2/3/4/5/7ml blood collection tube, 60*15ml(●), 60*15ml(▲), 24*50ml(▲), 24*50ml(●), 28*50ml(●), 28*50ml(▲)

2. Speed and centrifugal force

Compare the speed and relative centrifugal force your protocol requires with the rotor maximums.

- **Max. Speed(rpm):** 4000
- **Max. RCF(xg):** 3399

3. Centrifuge compatibility

Verify that your floor-standing low-speed centrifuge accepts this rotor type before buying it.

- **Rotors Type:** Swing Rotor



Importadora Gyb Limitada
Padre orellana 1476
Santiago CL13101 Chile

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)