

Rotor angular 12 x 20 ml para centrifuga de piso de baja velocidad

Código: B2-01-012213-R30

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Angle Rotor

Especificaciones técnicas

Especificación	-R30
Tipo de rotor	Rotor angular
Volumen	12*20ml
Velocidad máxima (rpm)	5000
Fuerza centrífuga relativa máxima (xg)	3516
Adaptadores	15ml, 10ml, 5ml, 1.5/2ml

Aplicaciones

1. Capacity and tubes

Confirm how many samples you process per run and the volume of each tube, and compare it with the rotor capacity.

- **Volume:** 12*20ml

Check whether your smaller tubes need one of the listed adapters.

- **Adapters:** 15ml, 10ml, 5ml, 1.5/2ml

2. Speed and centrifugal force

Compare the speed and relative centrifugal force your protocol requires with the rotor maximums.

- **Max. Speed(rpm):** 5000
- **Max. RCF(xg):** 3516

3. Centrifuge compatibility

Verify that your floor-standing low-speed centrifuge accepts this rotor type before buying it.

- **Rotors Type:** Angle Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)