

Rotor angular 12 x 15 ml para centrifuga de piso de baja velocidad

Código: B2-01-012213-R31

Marca: N/D

[Ver producto en Cotecno](#) | [Solicitar cotización](#)

[Detalle](#) | [Especificaciones](#) | [Aplicaciones](#)

Detalle

Angle Rotor

Especificaciones técnicas

Especificación	-R31
Tipo de rotor	Rotor angular
Volumen	12*15ml
Velocidad máxima (rpm)	5000
Fuerza centrífuga relativa máxima (xg)	3354
Adaptadores	1.5/2ml, tubo de recolección de sangre 2/3/4ml

Aplicaciones

1. Confirm tube capacity

Compare the number and volume of tubes you process per run with the rotor capacity.

- **Volume:** 12*15ml

2. Check speed and centrifugal force

Check that your protocol does not require a speed or relative centrifugal force above the rotor maximums.

- **Max. Speed(rpm):** 5000
- **Max. RCF(xg):** 3354

3. Check adapters and centrifuge

Confirm that the tubes you use match the adapters listed for this rotor.

- **Adapters:** 1.5/2ml, 2/3/4ml blood collection tube

Check that your centrifuge accepts this rotor type before purchase.

- **Rotors Type:** Angle Rotor

Siguiente paso: [ver producto en Cotecno](#) | [solicitar cotización](#)