

**PROBETA DE VIDRIO BELLTRONIC DE 250 ML CON BASE HEXAGONAL DE VIDRIO,
BELLPROBVID-250ML**



SKU: B-01-47-0106-1005 | **Categorías:** [Insumos de Laboratorio](#), [Material de vidrio de laboratorio](#), [Probeta de vidrio](#) |

DESCRIPCIÓN DEL PRODUCTO

| Codigo | Código SKU | Capacidad (ml) | Sub. División (ml) | Tolerancia + ml | Piezas/Contendor |
|--------------------|------------------------|----------------|--------------------|-----------------|------------------|
| BellProbVid-5ml | B-01-47-0106-1000 5 | | 0.1 | 0.1 | 360 |
| BellProbVid-10ml | B-01-47-0106-1001 10 | | 0.2 | 0.2 | 288 |
| BellProbVid-25ml | B-01-47-0106-1002 25 | | 0.5 | 0.5 | 288 |
| BellProbVid-50ml | B-01-47-0106-1003 50 | | 1 | 1 | 144 |
| BellProbVid-100ml | B-01-47-0106-1004 100 | | 1 | 1 | 96 |
| BellProbVid-250ml | B-01-47-0106-1005 250 | | 2 | 2 | 48 |
| BellProbVid-500ml | B-01-47-0106-1006 500 | | 5 | 5 | 24 |
| BellProbVid-1000ml | B-01-47-0106-1007 1000 | | 10 | 10 | 12 |
| CWGPV00001000 | CWGPV00001000 1000 | | - | - | - |
| BellProbVid-2000ml | B-01-47-0106-1008 2000 | | 20 | 20 | 8 |

INFORMACIÓN ADICIONAL

COTECNO