

CONJUNTO DE PRESIÓN DE NITRÓGENO A ALTA PRESIÓN BELLTRONIC



SKU: B-01-16-02-14-0500 | **Categorías:** [Filtro Prensa \(HTHP, Dinámico\)](#) |

DESCRIPCIÓN DEL PRODUCTO

[vc_row type="in_container" full_screen_row_position="middle" column_margin="default" scene_position="center" text_color="dark" text_align="left" overlay_strength="0.3" shape_divider_position="bottom" bg_image_animation="none"] [vc_column column_padding="no-extra-padding" column_padding_position="all" background_color_opacity="1" background_hover_color_opacity="1" column_link_target="_self" column_shadow="none" column_border_radius="none" width="1/1" tablet_width_inherit="default" tablet_text_alignment="default" phone_text_alignment="default" overlay_strength="0.3" column_border_width="none" column_border_style="solid" bg_image_animation="none"] [vc_column_text]**Conjunto de Presión de Nitrógeno a Alta Presión Belltronic B-01-16-02-14-0500**

[/vc_column_text] [/vc_column] [/vc_row] [vc_row type="in_container" full_screen_row_position="middle" column_margin="default" scene_position="center" text_color="dark" text_align="left" overlay_strength="0.3" shape_divider_position="bottom" bg_image_animation="none"] [vc_column column_padding="no-extra-padding" column_padding_position="all" background_color_opacity="1" background_hover_color_opacity="1" column_link_target="_self" column_shadow="none" column_border_radius="none" width="1/1" tablet_width_inherit="default" tablet_text_alignment="default" phone_text_alignment="default" overlay_strength="0.3" column_border_width="none" column_border_style="solid" bg_image_animation="none"] [vc_column_text]**Descripción del Conjunto de Presión de Nitrógeno a Alta Presión Belltronic B-01-16-02-14-0500:**

Conjunto de presión de nitrógeno a alta presión

Capacidad: 2 L[/vc_column_text] [/vc_column] [/vc_row]

INFORMACIÓN ADICIONAL

Volumen (ml)

[2000](#)

COTECNO