

EATON LOFMEM WB LMWB-10-1.2-4S FILTER CARTRIDGE



SKU: LMWB-10-1.2-4S DuP | **Categorías:** [Equipo Analítico](#) | **Etiquetas:** [brand|Lenntech](#)

DESCRIPCIÓN DEL PRODUCTO

EATON LOFMEM WB LMWB-10-1.2-4S FILTER CARTRIDGE

PART NUMBER: LMWB-10-1.2-4S

[Please contact us for more information and/or a quotation request.](#)

Specifications

Brand	Eaton LOFMEM WB
Model	LMWB 10 1.2 4S
Pore Size μ	1.2

Dimensions / Weights

Length in inch	10.0
Length in cm	25.4
Quantity per case	12 pieces



[Download Datasheet](#)

INFORMACIÓN ADICIONAL

url_fab	https://www.lenntech.com/products/Eaton-LOFMEM-WB/LMWB-10-1.2-4S/LMWB-10-1.2-4S-Filter-Cartridge/index.html
brand	Lenntech

COTECNO