

PORTAELECTRODOS - MARCA BELLTRONIC



SKU: N / A | **Categorías:** [Portaelectrodos](#) |

VARIACIONES

Imagen	SKU	Descripción
 A simple laboratory stand with a vertical rod and a horizontal arm for holding an electrode.	BellProBasE-600	Soporte básico de electrodos
 A laboratory stand with a circular base and a vertical rod, designed to hold an electrode and a stirrer.	BellProBasE-601	Portaelectrodo con Agitador
 A laboratory stand with a circular base and a flexible, articulated arm for holding an electrode.	BellProBasE-602	Portaelectrodos flexible

INFORMACIÓN ADICIONAL

COTECNO