

**PROBETA MEDIDORA DE PLÁSTICO, BASE HEXAGONAL 1000ML, BELLPLA-  
PROBINAMPBH-1000**



**Cod Alt: BellProbPI-1000**

**SKU:** B-01-21-0111-0404 | **Categorías:** [probetas de plastico](#) |

## DESCRIPCIÓN DEL PRODUCTO

| Modelo                   | Código Sku        | Tamaño disponible ml | Material      |
|--------------------------|-------------------|----------------------|---------------|
| BellPla-ProbinaMPBH-50   | B-01-47-0111-0400 | 50                   | polipropileno |
|                          | B-01-21-0111-0400 |                      |               |
| BellPla-ProbinaMPBH-100  | B-01-47-0111-0401 | 100                  | polipropileno |
|                          | B-01-21-0111-0401 |                      |               |
| BellPla-ProbinaMPBH-250  | B-01-47-0111-0402 | 250                  | polipropileno |
|                          | B-01-21-0111-0402 |                      |               |
| BellPla-ProbinaMPBH-500  | B-01-47-0111-0403 | 500                  | polipropileno |
|                          | B-01-21-0111-0403 |                      |               |
| BellPla-ProbinaMPBH-1000 | B-01-21-0111-0404 | 1000                 | polipropileno |
| BellPla-ProbinaMPBH-2000 | B-01-21-0111-0405 | 2000                 | polipropileno |

## INFORMACIÓN ADICIONAL

COTECNO