DISPOSABLE MICROBIOLOGICAL FILTER FUNNELS

These microbiological monitors are low-cost, disposable, presterilized units capable of both filtering and incubating samples while eliminating preparation and assembly time. The “all-in-one” sterile construction of these filter funnels makes them ideal for microbiological analysis tests, testing of food/beverages, raw material tests, sterility testing, and examining finished product quality. Because no flaming is required, and the filtration unit converts easily into a petri dish, there is minimal risk of crosscontamination.



FEATURES

- Pre-sterilized
- Cellulose nitrate membrane
- 0.45 µm pore size
- 47 mm diameter; 100 mL capacity

 Copyright © 2025