
Portable Suction System

Perfect size
Suitable for in-hood operations.

Storage rack
Keep the bench top tidy and spacious by accommodating all suction tools in storage.


PP vacuum bottle, 1000 mL
Optional 500 mL glass bottles are also available.
(Catalog number: 197000-60-500)



Liquid collection tray
Collects liquid waste dripping from suction tools and can be cleaned easily.

Fence-like platform
Tripping-proof and spill-proof.

Ergonomic handle
Lock switch keeps continuous suction and relives your finger muscle from repetitive procedures.



More energy-saving & quiter
Eco modles pause and resume operation automatically.

One-finger tip ejection
Effortless tip removal.

The product can be disinfected with 70% Ethanol or isopropanol (IPA).